

Portable Medical Ambulance Outside Ventilator/Handheld

Ambulance Ventilator YJ-PV150



Portable Medical Ambulance Outside Ventilator/Handheld Ambulance Ventilator

The Portable Medical Ambulance Outside Ventilator/Handheld Ambulance Ventilator is used for First-Aid, Ambulance, Emergency situation and patients' transport in hospital. It can be used for child and adult patients.

Product introduction of medical ventilator YJ-PV150:

Multiple functions, compact design, easy to take, design for transport and first aid. The main parts adopt good-quality components, accuracy and reliable. Special self-test function. Large LCD screen, simple and intuitive. 3 kinds of power sources can be chosen: AC, DC and Internal Battery.

Product parameter(specification)of the medical ventilator YJ-PV150:

Modes	CMV, A/C, SIMV, CPAP
-------	----------------------

Tidal Volume	50~2000mL
Respiratory Rate	1~100bpm
FiO2	40%~100%
Trigger	-20~19hPa
Peak pressure	0~60 hPa
Max Minute Volume	≥22L/min
PEEP	0, 5~20hpa Adjustable
PIP	0~60hPa
Over Pressure Relief	≤ 72hPa
Sigh	0~10 per 100 breathes
Package Dimension	530*530*280mm
Weight	13kg (Including all the accessories)

Product feature and application of medical ventilator YJ-PV150:

1. Microcomputer controlled, constant volume, pressure triggered.
2. LCD monitor, multiparameter digital display.
3. Light and handy, mainframe weights 3 kg, the overall weight is less than 7 kg (with the oxygen cylinder and all the accessories).
4. AC, DC and build-in Li-battery (can work for more than 8 hours) are all available.
5. Can be used for both invasive and noninvasive, both children and adults.

Product detail of medical ventilator YJ-PV150:

Technical parameters:

- Ventilation mode: CMV, A/C, SIMV, Sigh, SPONT, etc.
- Tidal Volume: 50 ~ 1200ml, continuously adjustable, digital display.
- Respiratory Rate: 1 ~ 60bpm, continuously adjustable, digital display.
- Oxygen Concentration: 40 ~ 100%, continuously adjustable
- Synchronous trigger sensitivity: -10 ~ 0hpa
- PEEP: 0, 5 ~20hpa

Monitoring and Alarm Function:

- Multiple parameters monitoring and display functions.
- Multiple alarm and prompt functions.
- Multiple security protection functions.

Service Conditions:

- Power Supply: AC, DC, Build-in Li-battery.
- Gas Supply: Oxygen cylinder, Central oxygen supply.

Basic configuration:

Mainframe of ventilator, carrying bag, silicone breathing circuits, oxygen cylinder, pressure reducing valve, high-pressure oxygen tube, power adapter, etc.

Product qualification of the medical ventilator YJ-PV100:

CE, ISO

Delivery, Shipping, and Serving of medical ventilator YJ-PV150:

Shipment: We are shipping products after receive payment 7 work days

Delivery time: 7 work days

Package: Standard package

After sale services: We provide one year's warranty. We guarantee help customers to settle all of quality problem with free during warranty.

Latest news:

Henan Forever Co., LTD have 20 years trading experience in the Chinese medical equipment marketing. Now company decides to open up global market. The company bring latest medical equipment products to the world and try to meet various demand in the global medical marketing.

Related Products:



YJ-PV100

YJ-PV80



YJ-PV600

WARRANTY:

This equipment will be enjoyed 1 year warranty according to the international standard .We will responsible for repairing it or replace the related accessories free of charge if any malfunction happened for the equipment's own quality problem from the date of purchase. Please contact us if you need further information for this model !