

Specification

*Product Name	128 Elements Waterproof Dustproof Cow Fetal Reproductive Disease Wireless USB Ultrasound Scanner Veterinary Ultrasound Machine
*Application scenarios	Engineered for all aspects of bovine and equine reproductive imaging, as well as for equine tendon imaging.
*Elements	128
*Imaging Modes	B, M, Color, Power, PW
*Gray scale levels	256
*Optional Probes	Rectal linear probe /Rectal convex probe
*Transducer frequency	Rectal linear probe: 6.0-11.0MHz Rectal convex probe: 2.0-5.0MHz
*Digital Technology	Tissue harmonic Imaging Tissue Specific Imaging Panoramic focusing technology Frequency compounding Space compounding Wide-angle imaging Real-time Dynamic Aperture
*Viewing device	Tablet, Phone, Goggles
*Ultrasound platform	Android, IOS(reserved), Windows(reserved)
*Scanning depth	up to 310mm
*Cine loop	400 frames movie playback stored
*Waterproof level	IP67
*Connection Type	WIFI, USB Type-C
*Product Size	190x90x60 mm
*Weight	630g (With probe)

